

EchelonVAC

Building on the established hydraulic technology of Echelon, EchelonVAC incorporates the use of an elevated vacuum system, allowing for increased suspension and a healthier residual limb environment. The hydraulic motion of the ankle acts like a pump, creating a vacuum between the socket and residual limb. This makes for a more secure attachment of the prosthesis, reducing the amount of movement (pistoning) present and increasing the safety and comfort of the user.

This vacuum also creates a negative pressure around the residual limb. Numerous studies have shown this pressure difference to be extremely beneficial in the management of wounds such as ulcers, a common ailment of amputees; promoting healthier skin and tissue at the residuum. The vacuum also helps to reduce the amount of volume fluctuation that occurs in the limb, ensuring a correct socket fit is maintained and reducing the risk of pain.

Improvements in Clinical Outcomes using Echelon compared to ESR feet

Improvement in **SAFETY**

- Reduced risk of tripping and falls
 - Increased minimum toe clearance during swing phase^{1,2}
- Improving standing balance on a slope
 - 24-25% reduction in mean inter-limb centre-of-pressure root mean square (COP RMS)³

Improvement in **ENERGY CONSUMPTION**

- Reduced energy expenditure during walking
 - Mean 11.8% reduction in energy use on level ground, across all walking speeds⁴
 - Mean 20.2% reduction in energy use on slopes, across all gradients⁴
 - Mean 8.3% faster walking speed for the same amount of effort⁴

Improvement in **MOBILITY**

- Improved gait performance
 - Faster self-selected walking speed^{2,5-7}
 - Higher PLUS-M scores than FlexFoot and FlexWalk style feet⁸
- Improved ground compliance when walking on slopes
 - Increased plantarflexion peak during level walking, fast level walking and cambered walking⁹
 - Increased dorsiflexion peak during level walking, fast level walking and cambered walking⁹
- Less of a prosthetic “dead spot” during gait
 - Reduced aggregate negative COP displacement⁵
 - Centre-of-pressure passes anterior to the shank statistically significantly earlier in stance⁵
 - Increased minimum instantaneous COM velocity during prosthetic-limb single support phase⁵

- Reduced peak negative COP velocity⁷
- Reduced COP posterior travel distance⁷
- Improved ground compliance when walking on slopes
 - Increased plantarflexion range during slope descent¹⁰
 - Increased dorsiflexion range during slope ascent¹⁰

Improvement in **RESIDUAL LIMB HEALTH**

- Helps protect vulnerable residual limb tissue, reducing likelihood of damage
 - Reduced peak stresses on residual limb¹¹
 - Reduced stress RMS on residual limb¹¹
 - Reduced loading rates on residual limb¹¹

Improvement in **LOADING SYMMETRY**

- Greater contribution of prosthetic limb to support during walking
 - Increased residual knee negative work⁶
- Reduced reliance on sound limb for support during walking
 - Reduced intact limb peak hip flexion moment⁶
 - Reduced intact limb peak dorsiflexion moment⁶
 - Reduced intact ankle negative work and total work⁶
 - Reduced intact limb total joint work⁶
- Better symmetry of loading between prosthetic and sound limbs during standing on a slope
 - Degree of asymmetry closer to zero for 5/5 amputees³
- Reduced residual and sound joint moments during standing of a slope
 - Significant reductions in both prosthetic and sound support moments¹²
- Less pressure on the sole of the contralateral foot
 - Peak plantar-pressure¹³
- Improved gait symmetry
 - Reduced stance phase timing asymmetry¹⁴

Improvement in **USER SATISFACTION**

- Patient reported outcome measures indicate improvements
 - Mean improvement across all Prosthesis Evaluation Questionnaire domains¹⁵
 - Bilateral patients showed highest mean improvement in satisfaction¹⁵
- Subjective user preference for hydraulic ankle
 - 13/13 participants preferred hydraulic ankle¹³

Improvements in Clinical Outcomes using EVS compared to other suspension types

Improvement in **SAFETY**

- Fewer falls and less chance of multiple falls
 - No trans-tibial EVS users reported multiple falls, while 75% of the non-EV users did¹⁶
- Better balance in functional clinical tests

- Significant improvements in the Berg Balance Scale (BBS), the Four Square Step Test (FSST) and the Timed-Up-and-Go (TUG) test¹⁷
- Better balance reported in patient-reported outcome measures
 - Improvements in the Activity Balance Confidence (ABC) scale questionnaire¹⁸

Improvement in **MOBILITY**

- Fewer gait compensations¹⁹⁻²¹
- Knee contact forces not significantly different to those of able-bodied controls²²

Improvement in **SUSPENSION**

- Decreased pistoning
 - Reductions of over 69% and 83%, compared to suction^{21,23} and pin-lock²⁴ suspensions, respectively, with other researchers and practitioners reporting similar observations^{18,19,25,26}
- Maintain residual limb volume
 - Suction suspension = mean 6.5% loss in volume; EVS = mean 3.7% increase in volume (N.B. it is possible that the increase may have been due to the fact that these individuals attended the clinic wearing their regular prostheses before using the EVS system)²¹
 - Other studies have since confirmed the observation that residuum volume loss is prevented by EVS^{19,27-30}

Improvement in **RESIDUAL LIMB HEALTH**

- Healthier residual limb tissue and skin
 - Higher trans-cutaneous oxygen measurement after activity³¹
 - Decreased trans-epidermal water loss after activity³¹
 - Decreased attenuated reactive hyperemia³¹
- Reduced interface pressures
 - Pressures reduced by a mean of 4% compared to suction suspension³²
 - Pressure impulses reduced by a mean of 7.5% compared to suction suspension³²
- Improved wound management
 - Continued prosthesis use while the wounds healed³³⁻³⁵
 - Wounds heal more quickly with EVS than other suspension methods³⁶
- Less painful than other suspension methods
 - Expert opinion¹⁹ and clinical case studies³⁷ agree that EVS is less painful and more comfortable than other suspension methods.
 - Improved Socket Comfort Score compared to other suspension methods³⁸

Improvement in **USER SATISFACTION**

- Patients are more satisfied wearing their prosthesis^{18,19,26,34,37-38}

Other Internal unpublished Blatchford research

Vacuum levels generated:

When sensory control of the lower limb joints is lost, it is essential that the replacement behaves predictably. Consistency of performance is vital in providing prosthetic confidence. In terms of socket suspension method, this means providing the same good connection throughout a gait cycle, from one step to the next, and day-to-day, over the lifetime of the socket.

The difference between the vacuum levels generated by suction suspension, and that generated when using EVS, can be demonstrated by using a negative pressure gauge⁴⁰. Figure 1 illustrates these measurements. Commonly, when the user bears weight on their prosthesis during stance phase, with suction suspension, the magnitude of the vacuum is low. When the leg is lifted into swing phase, the vacuum increases in magnitude, holding the socket to the residual limb. Comparatively, EVS retains a high level during stance phase – higher, in fact, than the peak swing phase vacuum with suction. Additionally, the difference between stance and swing phase is less pronounced, so that the vacuum level is more consistent throughout the gait cycle. For the amputee illustrated in the graph⁴⁰, EVS gave an approximate 85% increase in peak vacuum magnitude and an approximate 67% reduction in the 'amplitude' of the vacuum measurement signal.




Figure 1: Negative pressure within the socket when walking using a one-way valve suction suspension (grey) and an elevated vacuum (EV) suspension. N.B. Data recorded with EchelonVAC system.

The difference in vacuum generated by the AvalonVAC, compared to that generated by the EchelonVAC, is shown in Figure 2. Despite differences in the method used (keel vs springs, different socket, different pressure gauge), when the same patient was asked to walk at 'K2 walking speed' (~2km/h, short steps), the trend of vacuum level to number of steps taken was comparable to when measured at 'K3 walking speed' (4-5km/h) with EchelonVAC.




Figure 2: Comparison of the EchelonVAC and AvalonVAC vacuum generation by number of steps (regardless of walking speed)

References

1. Riveras M, Ravera E, Ewins D, Shaheen AF, Catalfamo-Formento P. Minimum toe clearance and tripping probability in people with unilateral transtibial amputation walking on ramps with different prosthetic designs. *Gait & Posture*. 2020 Sep 1;81:41-8.
2. Johnson L, De Asha AR, Munjal R, et al. Toe clearance when walking in people with unilateral transtibial amputation: effects of passive hydraulic ankle. *J Rehabil Res Dev* 2014; 51: 429.
3. McGrath M, Laszczak P, Zahedi S, et al. Microprocessor knees with “standing support” and articulating, hydraulic ankles improve balance control and inter-limb loading during quiet standing. *J Rehabil Assist Technol Eng* 2018; 5: 2055668318795396.
4. Askew GN, McFarlane LA, Minetti AE, et al. Energy cost of ambulation in trans-tibial amputees using a dynamic-response foot with hydraulic versus rigid ‘ankle’: insights from body centre of mass dynamics. *J NeuroEngineering Rehabil* 2019; 16: 39.
5. De Asha AR, Munjal R, Kulkarni J, et al. Impact on the biomechanics of overground gait of using an ‘Echelon’ hydraulic ankle-foot device in unilateral trans-tibial and trans-femoral amputees. *Clin Biomech* 2014; 29: 728–734.
6. De Asha AR, Munjal R, Kulkarni J, et al. Walking speed related joint kinetic alterations in trans-tibial amputees: impact of hydraulic ‘ankle’ damping. *J Neuroengineering Rehabil* 2013; 10: 1.
7. De Asha AR, Johnson L, Munjal R, et al. Attenuation of centre-of-pressure trajectory fluctuations under the prosthetic foot when using an articulating hydraulic ankle attachment compared to fixed attachment. *Clin Biomech* 2013; 28: 218–224.
8. Wurdeman SR, Stevens PM, Campbell JH. Mobility analysis of Amputees (MAAT 5): Impact of five common prosthetic ankle-foot categories for individuals with diabetic/dysvascular amputation. *J Rehabil Assist Technol Eng* 2019; 6: 2055668318820784.
9. Bai X, Ewins D, Crocombe AD, et al. Kinematic and biomimetic assessment of a

- hydraulic ankle/foot in level ground and camber walking. PLOS ONE 2017; 12: e0180836.
10. Bai X, Ewins D, Crocombe AD, et al. A biomechanical assessment of hydraulic ankle-foot devices with and without micro-processor control during slope ambulation in trans-femoral amputees. PLOS ONE 2018; 13: e0205093.
 11. Portnoy S, Kristal A, Gefen A, et al. Outdoor dynamic subject-specific evaluation of internal stresses in the residual limb: hydraulic energy-stored prosthetic foot compared to conventional energy-stored prosthetic feet. Gait Posture 2012; 35: 121–125.
 12. McGrath M, Davies KC, Laszczak P, et al. The influence of hydraulic ankles and microprocessor-control on the biomechanics of trans-tibial amputees during quiet standing on a 5° slope. Can Prosthet Orthot J; 2.
 13. Moore R. Effect of a Prosthetic Foot with a Hydraulic Ankle Unit on the Contralateral Foot Peak Plantar Pressures in Individuals with Unilateral Amputation. JPO J Prosthet Orthot 2018; 30: 165–70.
 14. Moore R. Effect on Stance Phase Timing Asymmetry in Individuals with Amputation Using Hydraulic Ankle Units. JPO J Prosthet Orthot 2016; 28: 44–48.
 15. Sedki I, Moore R. Patient evaluation of the Echelon foot using the Seattle Prosthesis Evaluation Questionnaire. Prosthet Orthot Int 2013; 37: 250–254.
 16. Rosenblatt NJ, Ehrhardt T. The effect of vacuum assisted socket suspension on prospective, community-based falls by users of lower limb prostheses. Gait Posture, <http://www.sciencedirect.com/science/article/pii/S096663621730111X> (2017, accessed 2 May 2017).
 17. Samitier CB, Guirao L, Costea M, et al. The benefits of using a vacuum-assisted socket system to improve balance and gait in elderly transtibial amputees. Prosthet Orthot Int 2016; 40: 83–88.
 18. Ferraro C. Outcomes study of transtibial amputees using elevated vacuum suspension in comparison with pin suspension. JPO J Prosthet Orthot 2011; 23: 78–81.
 19. Gholizadeh H, Lemaire ED, Eshraghi A. The evidence-base for elevated vacuum in lower limb prosthetics: Literature review and professional feedback. Clin Biomech 2016; 37: 108–116.
 20. Xu H, Greenland K, Bloswick D, et al. Vacuum level effects on gait characteristics for unilateral transtibial amputees with elevated vacuum suspension. Clin Biomech Bristol Avon 2017; 43: 95–101.
 21. Board WJ, Street GM, Caspers C. A comparison of trans-tibial amputee suction and vacuum socket conditions. Prosthet Orthot Int 2001; 25: 202–209.
 22. Xu H, Greenland K, Bloswick D, et al. Vacuum level effects on knee contact force for unilateral transtibial amputees with elevated vacuum suspension. J Biomech 2017; 57: 110–116.
 23. Gerschutz MJ, Hayne ML, Colvin JM, et al. Dynamic Effectiveness Evaluation of

- Elevated Vacuum Suspension. *JPO J Prosthet Orthot* 2015; 27: 161–165.
24. Klute GK, Berge JS, Biggs W, et al. Vacuum-assisted socket suspension compared with pin suspension for lower extremity amputees: effect on fit, activity, and limb volume. *Arch Phys Med Rehabil* 2011; 92: 1570–1575.
 25. Darter BJ, Sinitski K, Wilken JM. Axial bone-socket displacement for persons with a traumatic transtibial amputation: The effect of elevated vacuum suspension at progressive body-weight loads. *Prosthet Orthot Int* 2016; 40: 552–557.
 26. Scott H, Hughes J. Investigating The Use Of Elevated Vacuum Suspension On The Adult PFFD Patient: A Case Study. *ACPOC* 2013; 19: 7–12.
 27. Youngblood RT, Brzostowski JT, Hafner BJ, et al. Effectiveness of elevated vacuum and suction prosthetic suspension systems in managing daily residual limb fluid volume change in people with transtibial amputation. *Prosthet Orthot Int* 2020; 0309364620909044.
 28. Sanders JE, Harrison DS, Myers TR, et al. Effects of elevated vacuum on in-socket residual limb fluid volume: Case study results using bioimpedance analysis. *J Rehabil Res Dev* 2011; 48: 1231.
 29. Street G. Vacuum suspension and its effects on the limb. *Orthopadie Tech* 2006; 4: 1–7.
 30. Goswami J, Lynn R, Street G, et al. Walking in a vacuum-assisted socket shifts the stump fluid balance. *Prosthet Orthot Int* 2003; 27: 107–113.
 31. Rink C, Wernke MM, Powell HM, et al. Elevated vacuum suspension preserves residual-limb skin health in people with lower-limb amputation: Randomized clinical trial. *J Rehabil Res Dev* 2016; 53: 1121–1132.
 32. Beil TL, Street GM, Covey SJ. Interface pressures during ambulation using suction and vacuum-assisted prosthetic sockets. *J Rehabil Res Dev* 2002; 39: 693.
 33. Hoskins RD, Sutton EE, Kinor D, et al. Using vacuum-assisted suspension to manage residual limb wounds in persons with transtibial amputation: a case series. *Prosthet Orthot Int* 2014; 38: 68–74.
 34. Traballesi M, Delussu AS, Fusco A, et al. Residual limb wounds or ulcers heal in transtibial amputees using an active suction socket system. A randomized controlled study. *Eur J Phys Rehabil Med* 2012; 48: 613–23.
 35. Traballesi M, Averna T, Delussu AS, et al. Trans-tibial prosthesis in large area of residual limb wound: Is it possible? A case report. *Disabil Rehabil Assist Technol* 2009; 4: 373–375.
 36. Brunelli S, Averna T, Delusso M, et al. Vacuum assisted socket system in transtibial amputees: Clinical report. *Orthop-Tech Q Engl Ed*; 2.
 37. Arndt B, Caldwell R, Fatone S. Use of a partial foot prosthesis with vacuum-assisted suspension: A case study. *JPO J Prosthet Orthot* 2011; 23: 82–88.
 38. Rosenblatt NJ, Ehrhardt T, Fergus R, et al. Effects of Vacuum-Assisted Socket Suspension on Energetic Costs of Walking, Functional Mobility, and Prosthesis-

Related Quality of Life. *JPO J Prosthet Orthot* 2017; 29: 65–72.

39. Carvalho JA, Mongon MD, Belangero WD, et al. A case series featuring extremely short below-knee stumps. *Prosthet Orthot Int* 2012; 36: 236–238.
40. Sutton E, Hoskins R, Fosnight T. Using elevated vacuum to improve functional outcomes: A case report. *JPO J Prosthet Orthot* 2011; 23: 184–189.
41. McGrath M, Laszczak P, McCarthy J, et al. The biomechanical effects on gait of elevated vacuum suspension compared to suction suspension. Cape Town, South Africa, 2017.