

Multiflex

Multiflex uses a flexible articulation design attached to a foot with a keel that allows plantarflexion and controlled tibial progression, all the while remaining durable. It is adjustable in increments of approximately 5mm throughout its range via serrations at the interface of the Multiflex foot and the ankle.

Multiflex Slim also uses a flexible articulation design attached to a foot with a keel that allows plantarflexion and controlled tibial progression, all the while remaining durable. The main difference with the standard Multiflex is the nominal heel height. It is higher than a standard Multiflex foot, so not suitable for lower heeled shoes. It is adjustable, like a standard Multiflex foot, in increments of approximately 5mm throughout its range via serrations at the interface of the Multiflex foot and the ankle.

Multiflex SE also uses a flexible articulation design attached to a foot with a keel that allows plantarflexion and controlled tibial progression, all the while remaining durable. It was designed to negate the need for an ankle fairing to be used, as is common with regular Multiflex feet. This allows a shin fairing to interface directly with the SE foot. It has the same ROM as a conventional Multiflex foot.

Multiflex Ankle Standard 30mm allows customisable planterflexion and dorsiflexion by the appropriate selection of the rubber elements within the ankle. It is designed for use in conjunction with Multiflex feet. Multiflex Pyramid Ankle has a male pyramid attachment, rather than a 30mm tube clamp adaptor.

Clinical Outcomes using Multiflex feet

Multiflex was the “habitual” foot for all or majority of participants in 13 different studies¹⁻¹³.

With respect to **SAFETY**

- Low stiffness at weight acceptance leads to early foot-flat and greater stability for lower mobility patients¹⁴
- No loss of stability during standing with Multiflex than fixed ankle/foot¹⁵
- Easier to walk on uneven ground with Multiflex than fixed ankle/foot^{15,16}
- Easier to walk up a slope with Multiflex than fixed ankle/foot¹⁵

With respect to **MOBILITY**

- Little to no difference in gait mechanics compared to flexible, “energy storing” prosthetic feet¹⁷
- Increased prosthetic ankle range-of-motion with Multiflex compared to fixed ankle/foot^{15,16,18-20}
- Increased prosthetic ankle power with Multiflex compared to fixed ankle/foot for bilateral users¹⁶
- Prosthetic rollover shape closer to natural biomechanics than fixed ankle/foot¹⁸
- Bilateral users can walk longer distances and report “smoother” gait with Multiflex compared to fixed ankle/foot¹⁶
- Benefits in mobility for bilateral users^{15,16,18,19}

With respect to **RESIDUAL LIMB HEALTH**

- Equivalent socket comfort to higher technology, energy-storing feet²¹

With respect to **LOADING SYMMETRY**

- Improved stance phase timing symmetry with Multiflex compared to fixed ankle/foot²⁰
- Reduced sound limb loading with Multiflex compared to fixed ankle/foot²⁰

With respect to **USER SATISFACTION**

- Majority of users rate Multiflex as either “good” or “acceptable”²² and bilateral users prefer Multiflex to fixed ankle/foot¹⁶

References

1. Moore R. Patient Evaluation of a Novel Prosthetic Foot with Hydraulic Ankle Aimed at Persons with Amputation with Lower Activity Levels. *JPO J Prosthet Orthot* 2017; 29: 44–47.
2. Moore R. Effect on Stance Phase Timing Asymmetry in Individuals with Amputation Using Hydraulic Ankle Units. *JPO J Prosthet Orthot* 2016; 28: 44–48.
3. Buckley JG, De Asha AR, Johnson L, et al. Understanding adaptive gait in lower-limb amputees: insights from multivariate analyses. *J Neuroengineering Rehabil* 2013; 10: 98.
4. Sedki I, Moore R. Patient evaluation of the Echelon foot using the Seattle Prosthesis Evaluation Questionnaire. *Prosthet Orthot Int* 2013; 37: 250–254.
5. Rogers JP, Strike SC, Wallace ES. The effect of prosthetic torsional stiffness on the golf swing kinematics of a left and a right-sided trans-tibial amputee. *Prosthet Orthot Int* 2004; 28: 121–131.
6. Kobayashi T, Orendurff MS, Boone DA. Dynamic alignment of transtibial prostheses through visualization of socket reaction moments. *Prosthet Orthot Int* 2015; 39: 512–516.
7. Wright D, Marks L, Payne R. A comparative study of the physiological costs of walking in ten bilateral amputees. *Prosthet Orthot Int* 2008; 32: 57–67.
8. Vanicek N, Strike SC, Polman R. Kinematic differences exist between transtibial amputee fallers and non-fallers during downwards step transitioning. *Prosthet Orthot Int* 2015; 39: 322–332.
9. Barnett C, Vanicek N, Polman R, et al. Kinematic gait adaptations in unilateral transtibial amputees during rehabilitation. *Prosthet Orthot Int* 2009; 33: 135–147.
10. Emmelot C, Spauwen P, Hol W, et al. Case study: Trans tibial amputation for reflex sympathetic dystrophy: Postoperative management. *Prosthet Orthot Int* 2000; 24: 79–82.
11. Boonstra A, Fidler V, Eisma W. Walking speed of normal subjects and amputees: aspects of validity of gait analysis. *Prosthet Orthot Int* 1993; 17: 78–82.
12. Datta D, Harris I, Heller B, et al. Gait, cost and time implications for changing from PTB

- to ICEX® sockets. *Prosthet Orthot Int* 2004; 28: 115–120.
13. de Castro MP, Soares D, Mendes E, et al. Center of pressure analysis during gait of elderly adult transfemoral amputees. *JPO J Prosthet Orthot* 2013; 25: 68–75.
 14. Major MJ, Scham J, Orendurff M. The effects of common footwear on stance-phase mechanical properties of the prosthetic foot-shoe system. *Prosthet Orthot Int* 2018; 42: 198–207.
 15. McNealy LL, A. Gard S. Effect of prosthetic ankle units on the gait of persons with bilateral trans-femoral amputations. *Prosthet Orthot Int* 2008; 32: 111–126.
 16. Su P-F, Gard SA, Lipschutz RD, et al. Gait characteristics of persons with bilateral transtibial amputations. *J Rehabil Res Dev* 2007; 44: 491–502.
 17. Boonstra A, Fidler V, Spits G, et al. Comparison of gait using a Multiflex foot versus a Quantum foot in knee disarticulation amputees. *Prosthet Orthot Int* 1993; 17: 90–94.
 18. Gard SA, Su P-F, Lipschutz RD, et al. The Effect of Prosthetic Ankle Units on Roll-Over Shape Characteristics During Walking in Persons with Bilateral Transtibial Amputations. *J Rehabil Res Dev* 2011; 48: 1037.
 19. Major MJ, Stine RL, Gard SA. The effects of walking speed and prosthetic ankle adapters on upper extremity dynamics and stability-related parameters in bilateral transtibial amputee gait. *Gait Posture* 2013; 38: 858–863.
 20. Van der Linden M, Solomonidis S, Spence W, et al. A methodology for studying the effects of various types of prosthetic feet on the biomechanics of trans-femoral amputee gait. *J Biomech* 1999; 32: 877–889.
 21. Graham LE, Datta D, Heller B, et al. A comparative study of conventional and energy-storing prosthetic feet in high-functioning transfemoral amputees. *Arch Phys Med Rehabil* 2007; 88: 801–806.
 22. Mizuno N, Aoyama T, Nakajima A, et al. Functional evaluation by gait analysis of various ankle-foot assemblies used by below-knee amputees. *Prosthet Orthot Int* 1992; 16: 174–182.