# Tectus Evaluation Brace

## Instructions for Use

SU-TEB-L, SU-TEB-R

EN	Instructions for Use	2
DE	Gebrauchsanweisung	34
NO	Brugsanvisning	66

Contents .....	2
Packing Case (Peli™ 1640) Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	7
3 Construction .....	9
3.1 Component Identification and Adjustable Features .....	9
3.2 Height Adjustment Ranges .....	12
3.3 Cuff and Carbon Cushion Width Adjustment Ranges.....	14
3.4 Carbon Cushions .....	14
4 Brace Fitting and Alignment .....	15
4.1 Device Joint Damage Prevention.....	15
4.2 Alignment—Comparison with Custom-Made Tectus.....	16
4.3 Alignment Objectives.....	17
4.4 Device Set-up and Alignment Procedure.....	17
5 Tectus Knee Module Calibration and Configuration .....	25
6 Device Operation .....	25
7 Resetting the Device for the Next User .....	25
8 Maintenance .....	26
9 Limitations on Use .....	27
10 Assembly Instructions (Remove and Replace Cuffs) .....	28
11 Technical Data .....	29
12 Ordering Information .....	30

## Warning

**Clinicians are required to have received training on Tectus KAFO and to have read and understood the Tectus Clinician’s Manual 938465 prior to setting up and operating the Tectus Evaluation Brace.**

938465

### Packing Case (Peli™ 1640) Contents

- |  |  |
|--|--|
| 1. Fully assembled Tectus Evaluation Brace (including Tectus knee module and battery pack) | 7. Set of three Tectus knee module compression springs (stiffest one is pre-installed in the device) |
| 2. Free brace foot   | 8. Key fobs & lanyards (2 sets)  |
| 3. Packing cushions (option)   | 9. Hex key   |
| 4. Instructions for Use  | 10. Torx screwdriver   |
| 5. Third party ankle manual  | 11. 4-pin installation tool  |
| 6. Battery charger   | 12. Set of large lower thigh and calf cuffs  |

---

# 1 Description and Intended Purpose

## Preliminary Information

These instructions for use are for the practitioner who has completed Tectus training and is aware of the selection criteria for the patient. There is no separate IFU (Instructions for Use) document for the user (patient/wearer).

Ensure that you (the practitioner) understand all instructions for use, drawing particular attention to the safety information in Section 2 and cleaning and disinfecting instructions in Section 8. Ensure you have also read and understood the Tectus Clinician's Manual 938465, prior to using this device.

Inform the user of all the safety information relevant to them in Section 2 of this IFU, and instruct them that they must observe and comply with the safety warnings and advice given throughout the evaluation (i.e. while using the device and during fitting and removal).

## Definition of Terms

The term *device* is used throughout these instructions for use to refer to Tectus Evaluation Brace, an evaluation tool used in determining whether or not a patient may be suitable for being prescribed a custom-made Tectus.

The term *programming device* is used throughout these instructions for use to refer to a tablet or mobile phone running an Android or iOS operating system.

The term *app* is used throughout these instructions for use to refer to Tectus Programming Application which resides on the programming device.

The device includes an adjustable brace that operates like a Tectus KAFO (Knee Ankle Foot Orthosis) which can be adapted to a wide range of user leg and feet dimensions, hereafter referred to as *brace*.

The brace (see Fig. 2) forms a frame, hinged at the knee and ankle, that fits around the user's affected leg, and includes:

- Adjustable height side-bars to adjust device knee height
- Adjustable (A/P (Anterior/Posterior) Shift, rotation and width) cuffs with straps that fit around the thigh and calf
- A foot with adjustable heel stop, adjustable straps that fit around the user's shoe, and adjustable width foot/ankle side-bars
- Tectus knee module providing lateral-side knee joint
- Tectus knee follower joint (providing medial-side knee joint)
- Tectus battery mounting, Tectus Battery pack and power cable
- Third party ankle unit (lateral-side) and ankle follower joint (medial-side)

In addition to the brace, the device includes :

- A separate brace shoe with adjustable straps for the foot on the unaffected leg
- Tectus battery charger
- Two Tectus key fobs each with a lanyard
- Tectus Programming App (downloadable to Android or iOS programming device)
- A set of large lower thigh cuffs and calf cuffs (the brace is fitted with a standard size set).
- Two alternative compression springs for the Tectus module and associated fitting tools

## Scope

This document describes primarily the operation and setup of the brace, not the operation of the Tectus knee module or its associated app (which provides the primary system set-up and control adjustment functions of the device).

For information regarding the set-up, operation and control of the Tectus module, battery and key fob, or programming app, please refer to the Tectus Clinician's Manual (938465) and User Guide (938464).

Maintenance of the device is the responsibility of Blatchford in line with the terms of the loan agreement.

- Notes:
- i) In the related Tectus Clinician's Manual, the term *device* refer to the Tectus knee module.
  - ii) The device's ankle joints are manufactured by a third party, please refer to ankle joint manufacturer's documentation for their Instructions for use.

## Application

**The device is solely for use as a tool for assessing whether or not a patient could benefit from being prescribed a custom-made Tectus KAFO (Knee Ankle Foot Orthosis); it is not a medical device.**

**The device is for temporary use only. It must only be used under clinical supervision, and always used with parallel bars for support. It must not be used outside of a clinically controlled environment.**

The device is for use by multiple users.

The device is an adjustable, microprocessor-controlled, swing and stance, KAFO-like tool for simulating the function and operation of a custom-made Tectus KAFO. The device has potentially the same operation and controls as a custom-made Tectus, but with fully-adjustable brace components that allow it to be adapted to fit a wide range of users.

Left and right leg brace variants are available. For leg length equilibrium, each brace is supplied with a separate free foot with adjustable sole, heel-stop and fitting straps for the unaffected limb which must be worn at the same time as the brace on the affected limb during the evaluation.

In brief, the device is used as follows: The practitioner adjusts the brace knee height, leg width and heel-stop position and the number of fitted shoes blocks to suit the dimensions of a specific user and their footwear, then fits the brace to the user with loosely fastened straps. They then carefully, correctly align the user's anatomical compromise knee center A/P position relative to the brace knee center A/P position including any calculated offset, by horizontally sliding and/or rotating the lower thigh cuff, (upper thigh cuff if fitted) and calf cuff to their appropriate positions and set any required tilt of the shin, before locking the brace components in position by tightening their adjustment screws, then securely fastening and tightening all brace straps around the user's leg and shoe. The practitioner also adjusts the free foot to match the brace foot strapping and foot block arrangement and fits this over the user's shoe on their unaffected leg.

The practitioner then calibrates the integrated Tectus module and adjusts its program parameters to tailor device operation to suit the needs of the user, via the Tectus programming app. For example, they set how much flexion resistance the knee module's hydraulic cylinder will automatically apply during the swing gait phase; how much it will apply during sitting mode; how much in stairs mode, etc. They can also set the specific thigh angle (extension) that the Tectus knee module must detect for the knee to release into swing.

Note... Note that extension assistance (extension bias) is spring-controlled and cannot be adjusted using the app. The practitioner must manually adjust this, as necessary, by tightening or loosening the spring adjustment screw, or by swapping the spring for one which has a different stiffness. (Each device includes 3 springs of differing stiffness, the stiffest is fitted in the Tectus module by default.) For further information, refer to the Tectus Clinician's Manual 938465.

The user then uses the correctly set-up device as directed by the practitioner. This allows the user to experience, in a safe and controlled environment, how a Tectus would work for them; how it feels to use such a device for activities such as walking, standing, sitting and stepping up on or down from a step, and feedback their observations to the practitioner. At the same time it allows the practitioner to observe and record the user's gait performance as they use the device. The practitioner can then decide whether or not it would be beneficial to proceed in prescribing a custom-made Tectus KAFO for the user.

During the evaluation, the practitioner can switch the device between different user operating modes for different activities, using the app. The practitioner can, or user may (if requested to do so during the evaluation), also switch user modes using control buttons on the device, or using control buttons on a key fob that has been paired with the device.

- Notes:
- i) The device is worn on the user's affected leg and has a brace foot that fits around the outside of the user's footwear (unlike the Tectus custom brace where its brace foot sole fits inside the user's footwear).
  - ii) If the user currently uses an orthosis that fits inside their shoe, and therefore may wear a shoe larger than normal size, it is advised they bring a pair of normal sized shoes in addition so that they could try the Evaluation brace with or without their existing orthosis.

## Features

### Brace

- Slipper toe for adjustable foot length to fit shoe sizes 5–11 (UK), 38–46 (EU), 6–12 (US)
- Adjustable knee height (and adjustable thigh height if upper thigh cuff is fitted)
- Sagittal A/P shift-adjustable heel stop position
- Knee center sagittal adjustment via adjustable (A/P shift and rotation) calf cuff and lower thigh cuff
- Adjustable width at upper (if fitted) and lower thigh cuff, calf cuffs and foot
- Foot with quick-release cam buckle (cargo) straps and hook and loop (Velcro®) fastening straps
- Pliable, floating, thigh and calf carbon cushions, strap-mounted on medial side, with:
  - QuickFit® 'CLICK®' quick-release, self-locating magnetic buckles
  - QuickFit® 'BOA®' ratchet straps with micro-adjustment dial
- Standard (fitted) and large (in packing case) sets of thigh and calf cuffs
- Supplied in sturdy, wheeled, Peli™ packing case for safe and easy transportation

## **Tectus Knee Module Features and Controls**

Refer to Tectus Clinician's Manual 938465

## **Third Part Ankle Features and Controls**

Refer to third party ankle documentation.

## **Indications**

The device is intended solely for use as an assessment tool to indicate whether or not a patient would benefit from being prescribed a custom-made Tectus KAFO (Knee Ankle Foot Orthosis); it is not a medical device.

The brace is intended for temporary, short period use only, typically no more than 30 minutes.

The device is only be used under clinical supervision and the brace should always be used with parallel bars support for the user.

The one-size-fits-all approach to the design means it must be accepted that there will be joint alignment compromises and user comfort may be not optimal.

For all medical indications associated with the Tectus Knee Module, please refer to the Tectus Clinician's Manual 938465 for details.

## **Contraindications**

The device is not suitable for users weighing more than 100 kg (220 lb) or less than 50 kg (110 lb).

The device is not suitable for users with a knee height (when wearing shoes to be used with the device) of greater than 54 cm (1 ft 9¼ in) or less than 42.5 cm (1 ft 4¾ in).

Not suitable without using an AFO (Ankle-Foot-Orthosis) for users with weak ankles that would normally use an AFO.

For all medical contraindications associated with the Tectus Knee Module, please refer to the Tectus Clinician's Manual 938465 for details.


The responsibility for whether or not this evaluation tool can be used by a specific patient ultimately lies with the Tectus-trained Practitioner.

## 2 Safety Information

-  This warning symbol highlights important safety information which must be followed carefully.
-  The device is intended solely for use as a tool for assessing whether or not a patient could benefit from being prescribed a custom-made Tectus KAFO (Knee Ankle Foot Orthosis); it is not a medical device.
-  The device is intended for temporary, short period use only, typically no more than 30 minutes, and must only be used within a clinical environment.
-  The device must only be used under clinical supervision and should always be used with parallel bars support
-  Assembly, programming, maintenance and repair of the device (including handling and replacement of the battery unit) must only be carried out by a suitably qualified practitioner or technician that has attended an approved training course.
-  Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.
-  Avoid strong magnetic fields, radioactive environments, sources of electrical interference, atmospheres containing liquids and/or powders.
-  The device must not be worn during any scanning procedures (e.g. X-ray, CT or MRI scan).
-  The device is not intended for use when immersed in water or as a shower orthosis. If the orthosis comes into contact with water remove and dry immediately.
-  Only footwear approved by the practitioner that fits securely within the brace should be used with the device.
-  Contact Blatchford or your service provider for advice on whether the device can be used by a user fitted with a pacemaker or any other electronic medical device.
-  Ensure no water enters the charger plug socket and the cover is only removed during charging.
-  Any changes in the performance or function of the device e.g. restricted movement, non-smooth motion or unusual noises should be reported to Blatchford/ your service provider.  
To prevent damage to user's clothing take care not to trap clothing between the moving parts of the device.
-  If the critical service indicator is activated and cannot be reset stop using the device immediately and contact Blatchford or your service provider to arrange a service.
-  Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
-  **WARNING:** Portable RF (radio frequency) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 in) to any part of the device, including cables specified by the manufacturer; this may otherwise result in the degradation of performance of the device.
-  Be aware of finger trap hazards at all times.
-  When disinfecting the device, read and observe the safety warnings associated with the disinfectant used. Wear eye protection and protective gloves, and carry out disinfection in a well-ventilated area and/or use an air respirator.

 The device is not self-supporting in **Free** mode. Ensure the user (or device if not being worn) is supported before switching to **Free** mode.

 After fitting and aligning the device, but before the user attempts to use the device, ensure all device adjustment screws are tight, all Velcro® straps and cam buckles are securely fastened, all cuff strap magnetic CLICK® fasteners are correctly inserted and fastened as shown in Fig. 1, and all BOA® ratchet strap micro-adjustment dials are pushed in and tightened (by rotating clockwise) so the device is securely fitted to the user.

Magnetic CLICK®  
fastener

BOA® Ratchet Strap  
Micro-adjustment Dial

Fig. 1

Advise the user to observe the following warnings and advice while they are using the device:

 Report any discomfort immediately to the practitioner. If you feel any numbness/loss of sensations, the device must be removed immediately and not reapplied until normal sensations have returned, and then only reapply with caution avoiding overtightening the straps.

 Be aware at all times of the bulkiness of the brace and free foot so that you do not trip over your own feet or bump into nearby objects.

 Be aware of finger trap hazards at all times.

 Always use with support e.g. parallel bars, hand rail or chair.

 Use a step-to-step gait when descending steps. Leg-over-leg step descent should only be undertaken after you have been trained and approved to do so by the practitioner.

 When starting to walk, always lead with the unaffected leg. If walking is initiated leading with the braced limb, the device may not recognize the movement and as a safety measure will switch to **Lock** mode.

 Confirm that **Stairs** mode is engaged before attempting to descend steps.

 Confirm that **Sit** mode is engaged before attempting to sit.

 Do not adjust or tamper with the setup of the device.

 After continuous use, the knee module may become hot to the touch.

 A tubular bandage or stocking must be worn over the affected leg if the clothing worn would otherwise allow direct skin contact with the brace.

 Do not carry loads while using the device during the evaluation.

 To prevent damage to clothing, take care not to trap clothing in the moving parts of the device.

 Avoid any strangulation hazard that the key fob lanyard may present to the user when worn, such as catching the cord on a fixed object while walking.

 Keep the key fob on a lanyard separate from any other key fob to prevent accidental mode switching.

 The device is not designed for walking backwards. Do not attempt to walk backwards as the device may not support you. In **Walk** mode, when backwards motion is detected, the device will lock in extension to provide standing support only.

# 3 Construction

## Principal Parts

Item	Material	Item	Material
Side-bars	Titanium	Carbon Cushions	Carbon fiber
Foot-bars	Stainless Steel	Feet	Thermoplastics
Cuffs	Nylon (PA12), Perlon	Strap fasteners	Steel
Screws	Stainless Steel	Straps	Polymers

### 3.1 Component Identification and Adjustable Features

**Left Leg, Front/Side View  
Medial Side (Inside)**

- Lower Thigh Carbon Cushion
- Lower Thigh Cuff
- Lower Thigh Cuff A/P Shift and Rotation Adjustment Screws
- Medial Thigh Side-bar (fixed length)
- Tectus Knee Follower Joint
- Calf Cuff A/P Shift and Rotation Adjustment Screws
- Extension Lock Screws
- Lower Extendable Knee Side-bar and Scale
- Third Party Ankle Follower Joint
- Medial Ankle Side-bar/foot-bar
- Mid Foot Strap
- Fore Foot Strap
- Sole
- Foot Blocks
- Slipper Toe

**Left Leg, Front/Side View  
Lateral Side (Outside)**

- Lateral Upper Thigh Cuff and Strap with Lateral Upper Extendable Side-bar (optional fit)
- Lower Thigh Cuff, Upper Strap
- Tectus Knee Module
- BOA® Ratchet Strap Micro-adjustment Dial (and CLICK Magnetic Fastener not shown, see Fig. 1)
- Lower Thigh Cuff, Lower Strap
- Tectus Knee Joint
- Calf Carbon Cushion
- Calf Cuff, Upper (Shin) Strap
- BOA® Ratchet Strap Micro-adjustment Dial (and CLICK Magnetic Fastener not shown, see Fig. 1)
- Calf Cuff, Lower (Shin) Strap
- Third Party Ankle
- Heel (Ankle) Foot Strap
- Top Foot Strap and Cam Buckle



Fig. 2

Note: All foot straps are hook and loop (Velcro®) fastening except top foot strap which uses a cam buckle fastening

(Optional Fit) Upper Thigh Cuff, Sagittal (A/P) Position and Tilt Adjustment Screw. (Sets anterior/posterior position and tilt angle.)

Fig. 3 c

**Rear/Side View Lateral Side (Outside)**

- Lower Thigh Cuff, Width Adjustment Screw
- Lower Thigh Cuff, Width Scale
- Lower Thigh Carbon Cushion
- Calf Carbon Cushion
- Calf Cuff, Width Adjustment Screw
- Calf Cuff, Width Scale
- Ankle Joint (Third Party)
- Cam Buckle


Fig. 3 a


Fig. 3


Fig. 3 b

- Side-bar/Foot-bar Width Adjustment Screws (Under Sole)
- Heel-Stop Sagittal (A/P) Position Adjustment Screw (Under Sole)

Note... Lower thigh cuff vertical height is fixed relative to the knee joint, only the upper thigh cuff (if fitted) and calf cuff have vertical height adjustment relative to the knee joint.

**Rear/Side View Medial Side (Inside)**

- Lower Thigh Cuff, Sagittal (A/P) Position and Tilt Adjustment Screws.


Fig. 3 d

- Calf Cuff Sagittal (A/P) Position and Tilt Adjustment Screws.

- Slide-in Foot Blocks
- Ankle Follower Joint

**Heel-Stop Sagittal (A/P) Position Sliding Adjustment**




Fig. 3 e

## Sagittal A/P Shift Range

- (A) Upper Thigh Cuff (lateral side only): Max.  $\pm 25$  mm
- (B) and (C), Lower Thigh Cuff: Max.  $\pm 32.5$  mm
- Lower Calf Cuff (not shown): Max.  $\pm 32.5$  mm


Each cuff can be shifted anterior or posterior in the sagittal plane (guided by angle adjusters that can slide along slots in the cuffs), and/or be rotated when the associated A/P adjustment screws are loosened.

## Sagittal Rotation Range

- Upper Thigh Cuff (not shown) (lateral side): Max.  $\pm 15^\circ$
- (D) Lower Thigh Cuff: Max.  $\pm 7^\circ$
- Calf Cuff (not shown): Max.  $\pm 6^\circ$


## Free Foot (for unaffected leg)

Identical to brace foot but with no ankle side bars fitted.

Removable Foot Blocks

Foot blocks slide in from open side.

After setting the heel stop A/P position on brace foot (see Section 4.4.8) and adding or removing foot blocks, set the free foot heel stop A/P position and foot block arrangement to match before fitting the free foot on the user.

When fitting the brace foot (or free foot), firstly fasten the heel (ankle) strap followed by the mid foot and fore foot straps around the user's shoe. Finally, thread the top strap through the heel strap (or mid foot strap as necessary) and back through the cam buckle before pulling tight to bend the flexible slipper toe around the user's shoe, then fasten the cam buckle for a secure fit.


Note... To adjust cam buckle straps, press cam buckle lever to release strap, release lever to lock.

## Heel-Stop Adjustment Range (Brace foot and free foot)


- (E) Max.  $\pm 10$  mm

Note... To gain access to adjust the lateral side, lower thigh cuff, lower A/P adjustment screw, or to remove/replace battery:

1. Loosen the Tectus module clamp screw and rotate the Tectus module
2. Loosen as necessary the lower A/P adjustment screw using the hex key supplied.
3. On completion of lower thigh cuff adjustments or battery replacement, reposition Tectus module and tighten clamp screw.

## 3.2 Height Adjustment Ranges

Key	Description	Height/cm
A	Device ankle height above top of sole	12
B	Minimum device knee center height above top of sole	42.5
C	Maximum device knee center height above top of sole	54
D	Minimum medial upper thigh cuff height above top of sole	66
E	Maximum medial upper thigh cuff height above top of sole	97
F	Height from ground to top of sole	2
G	Maximum height adjustment of device knee center	11.5
H	Maximum height adjustment of top of medial upper thigh cuff (Optional fit)	10.5

**Maximum Extension at  
Calf and Thigh**


Fig. 7

Key	Description	Height/cm
J	Overall maximum height	96
K	Lower side-bar height adjustment: maximum vertical distance between device knee and ankle centers	44
L	Lower side-bar height adjustment: minimum vertical distance between device knee and ankle centers	32
M	Lateral upper side-bar height adjustment: maximum height of top of medial cuff above device knee center (Optional fit)	40
N	Lateral upper side-bar height adjustment: minimum height of top of medial cuff above device knee center (Optional fit)	31
O	Height above knee center to top of lower thigh carbon cuff	21

### 3.3 Cuff and Carbon Cushion Width Adjustment Ranges

#### Cuff Selection Table

Cuff	Size	Internal Width Including Cushion	
		Minimum	Maximum
Upper Thigh Cuff (Optional fit)	Standard	n/a	n/a
	Large	n/a	n/a
Lower Thigh Cuff	Standard	11.5 cm (4 1/2 in)	16.5 cm (6 1/2 in)
	Large	15.5 cm (6 1/16 in)	20.5 cm (8 1/16 in)
Calf Cuff	Standard	10.6 cm (4 3/16 in)	15.6 cm (6 1/8 in)
	Large	14.6 cm (5 3/4 in)	18.6 cm (7 3/16 in)

### 3.4 Carbon Cushions

#### Upper Thigh Carbon Cushion Optional Fit for users typically taller than 160 cm (5 ft 3 in).

The lateral upper thigh carbon cushion is attached to the lateral upper-side bar (if fitted) and is height adjustable. The upper thigh strap is threaded around the lateral upper thigh carbon cushion and is free floating on the medial side. The medial strap, if used, should be wrapped around thigh and positioned 5 cm below the groin (or as close as is comfortable) and secured with the hook and loop fastening.


#### Lower Thigh and Calf Carbon Cushions

The carbon cushions are suspended inside the cuffs, on the medial side only, by the straps threading through them.

The carbon cushions are flexible and designed to fit around a range of thigh/calf sizes. They are fitted with QuickFit® straps with self-locating magnetic CLICK® fasteners and BOA® micro-adjustable ratchets. These straps wrap the carbon cushions tightly around the limb.

To engage a BOA® strap ratchet, push in the BOA® dial, then to tighten the strap, rotate the dial clockwise. (The dial will only turn clockwise). The dial clicks as it turns, incrementally tightening the strap by 1 mm with each click. To release the BOA® strap ratchet, pull the BOA® dial outwards.

Any large gap between the carbon cushions and the sides of the calf cuff, should be filled with soft padding pieces.

# 4 Brace Fitting and Alignment

## 4.1 Device Joint Damage Prevention.

To prevent damage to the device joints, after any brace adjustment the medial and lateral sides of the device side-bars must remain parallel from top to bottom, not twisted or tapering, i.e. with the medial and lateral knee joint pivot axes coaxially aligned and horizontal, and with the medial and lateral ankle joint axes coaxially aligned and horizontal. The width (separation) is set to the user's upper thigh width measurement .

This means that any brace knee or ankle height adjustment made on the medial side must be replicated exactly on the lateral side, and any brace width adjustment made for knee width at the thigh, must be matched by making the same width adjustment at the calf and at the foot, thus keeping the sides of the brace parallel. Similarly any cuff A/P shift or rotation must be made equally on both sides.

**Warning, incorrect alignment may result in damage, wear or malfunction of the device.**

Upper and lower medial side-bars, Tectus follower joint and ankle follower joint vertically aligned.

Upper and lower lateral side-bars and ankle joint vertically aligned with Tectus knee joint vertical axis, parallel with ankle joint and side bar vertical axis.

### Ideal Knee Alignment

User's Anatomical Compromise Knee Center Horizontal Axis

Device Knee Joints' Pivot Center Horizontal Axis

Ankle Joints' Horizontal Axis



Fig. 9 Front View

Build Line

Lower Thigh Cuff A/P shift and rotation as required to replicate as near as possible Tectus bench alignment

Calf Cuff A/P shift and rotation as required to replicate as near as possible actual Tectus bench alignment

Fig. 10 Medial Side View

## 4.2 Alignment—Comparison with Custom-Made Tectus

### Tectus Evaluation Brace

Left leg view on medial side


Fig. 11

### Custom-Made Tectus

Right leg view on lateral side


Fig. 12

Note... For alignment in the sagittal (A/P) plane, always align the user's anatomical compromise knee center as close as is practicable to the optimum calculated position, anterior and superior to the device knee center position.

Knee alignment must always take priority over ankle alignment.

After aligning knees and ankles, it may then not be possible to set the ideal angle of tilt of the shin that can be achieved in a custom-made device.

Set anterior tilt of knee as required (up to maximum 10°).

- D** Horizontal A/P depth of the knee.
- Device mechanical center of rotation at knee/ankle.
- Anatomical compromise knee center/ankle center.

## 4.3 Alignment Objectives

Aim to replicate as close as possible the bench alignment of the Custom-Made Tectus Brace, see Tectus Clinician's Manual 938465 and as shown in Section 4.2. Noting the limitations of the one-size-fits-all approach to the device design, a geometrical perfect fit may not be possible. Some compromises must be expected and it must be accepted that occasionally there may be some degree of discomfort for the user for the short period of expected usage.

Bearing in mind the need to keep the brace sides parallel, set the brace width (separation) to the user's upper thigh width measurement. If larger cuffs are required exchange the standard cuff set for the large cuff set.

Alignment should be medially biased. 7–12° of toe/forefoot extension is recommended, which can be achieved by tightening the toe strap.

Add sufficient foot-blocks to cover the user's front foot shoe sole. The top foot strap needs to be threaded through the ankle strap or mid sole strap then returned through the cam buckle and pulled tight to bend the slipper toe around the front of the user's shoe for security.

Note that the device ankle joint and ankle follower joint center height and A/P position are not adjustable, however the device's heel stop adjustment may provide some A/P adjustment of the user's shoe to help align their ankle center (medial /lateral average) closer to the device ankle center.

The device's medial knee joint's center should be aligned as the priority, ideally as close as possible to the description in the Tectus Clinician's Manual 938465. Aligning the user and device ankles and setting the shin/calf tilt is less critical for the device to work for the short duration of the evaluation.

The alignment need not be optimal for the Evaluation Brace device to provide adequate indication of whether or not the user could benefit from being prescribed a custom-made brace . However, the practitioner must ensure that the alignment is safe for the user.

Note... For the Tectus Evaluation Brace, the user's shoe fits inside the brace shoe i.e. it sits on top of the brace sole plate, whereas for the Tectus Custom Brace the sole plate fits inside the user's shoe. This means the thickness of the user's shoe and heel height and any toe/forefoot extension provided by the shoe must be taken into account when following the Bench Alignment procedure in the Tectus Clinician's Manual 938465 and applying them to the Tectus Evaluation Brace.

There are left (SU-TEB-L) and right (SU-TEB-R) device leg variants of the device. Ensure you have the correct brace variant for the affected leg of the user before starting set-up and alignment.


## 4.4 Device Set-up and Alignment Procedure

Note... Metric tape measure/ruler, internal and external calipers, laser line and weighing scales are recommended to aid measurement and alignment.

The following alignment procedure is a suggested approach. Practitioners, after training, may chose a different approach depending on their training, available equipment and knowledge/ understanding of their patients' conditions.

For most of the fitting, it is recommended that the user does not bear weight through the device; however, while final alignment adjustments are being made to the fitting, the user will need to be able to stand for a short time while being supported e.g. holding on to a chair, rail or parallel bars, until this is completed. Experienced practitioners who are aware of the user's ability to stand may chose to make certain measurements or carry out alignment with the user standing supported, or sitting on an examination table or chair as appropriate.

## Device Set-up and Alignment Flow

## 4.4.1 Power On and Reset the Tectus Knee Module

1. Ensure the Tectus is fully charged.
2. Power on the Tectus Knee Module and connect the app (see Tectus Clinician's Manual).
3. Check the basic device settings as shown in the app's BASIC SETUP screen are:
  - Sitting Resistance (50%)
  - Pressure Release Threshold (20 Bar)
  - Stance Release Angle ( $-10^{\circ}$ )
  - Stairs Resistance (50%)
  - Swing Adjustment (50%):and reset as necessary to these default settings:
4. Fully extended the device and set to **Lock** mode.

Note... **Free** mode may be selected temporarily to facilitate donning and/or measurement.

**In Free mode the device will not self-support, it will be free to swing at the knee, so ensure that the user, when wearing the device, and also the device itself when not worn, are safely supported before switching to Free mode. Be aware of the potential finger trap hazards at all times, and in particular when the device is free to swing at the knee.**

## 4.4.2 Weigh and Measure the User

### Weigh the User

Check the weight of the user. They must weigh less than 100 kg and greater than 50 kg otherwise they cannot use the device and you cannot proceed. Measure the following user dimensions if these have not already been provided.

### Measure the User

1. Position the user as shown in Fig. 13, sitting comfortably on a horizontal examination table.
2. Ensure the user is wearing appropriate footwear and suitable clothing that they will wear when using the device.

**Note... A tubular bandage or stocking must be worn over the affected leg if the clothing worn would otherwise allow direct skin contact with the brace.**

3. Measure the following dimensions (refer to Fig. 13) with the affected leg lying horizontal, and include clothing thickness in the measurements:
  - A User's thigh width (medial-lateral) at 21 cm above their knee center (equivalent to the height above device knee center of the top edge of the lower thigh carbon cushion).
  - B User's horizontal (when standing) A-P depth at the knee.
  - C User's knee vertical height (when standing): Measure from the bottom face of shoe heel horizontally to tibial plateau.
  - D User inside leg measurement: Measure from the bottom face of shoe heel horizontally to 5 cm below groin.


Fig. 13

- Notes: i) The thigh width A, measured at 21 cm (8¼ in) vertically above knee center (equivalent to top of lower thigh cuff carbon cushion on the device) is needed to set the width of the lower thigh cuff, and calf cuff using the width scales on the cuffs. At the same time, the ankle/foot side-bars separation is adjusted to match, such that the sides of the brace remain parallel. The user's measured thigh width will be used to select which cuff size, standard or large, to use on the brace.
- ii) Measurements B (knee depth) and C (knee height to tibial plateau) are used to calculate the location of the user's anatomical compromise knee center, and then to determine the ideal position of device knee joint center relative to the position of the user's anatomical compromise knee center.
4. Check the measured dimensions are within the range of adjustment limits of the device. For height dimensions, see Section 3.2; for widths, check required cuff size using the *Cuff Selection Table* in Section 3.3. If larger cuff sets are required fit them as per Section 10. If the user dimensions are outside of these range limits, do not attempt to fit the device.

### Calculate Required Device Knee Center Offset

1. Calculate the vertical and horizontal offset required of the device knee and follower joint centers relative to the user's knee joint anatomical compromise knee center (see training notes).

#### 4.4.3 Set the Lower Thigh and Calf Cuff Widths

1. Slacken the medial and lateral foot-bar securing screws so they are free to slide open/close when the knee and thigh cuffs are adjusted.
2. Slacken all the cuff/carbon cushion straps.
3. Slacken the feet and ankle straps.
4. Loosen the two foot bar adjustment screws under the sole.
5. Loosen the lower thigh cuff and calf cuff width adjustment screws and slide the cuff sides together or apart equally as required to set the cuff width to the user's maximum thigh width A (see Fig. 13) using the centimeter scale on each cuff, then tighten each adjustment screw to 5 Nm torque to lock in position. See Fig. 14 for how to read the cuff scale.

Note... Cuff scales are read at the arrow position shown in Fig. 14, where the inner vertical edge of the outer cuff (directly below hole center) lines up with the scale on the inner cuff. Example shown is a lower thigh cuff width setting of 14.5 cm.

**It is critical that the thigh and calf cuff widths are set to the same width for the device to function correctly and to prevent damage to the knee and ankle joints.**


Fig. 14

#### 4.4.4 Set the Foot Bar Separation to Maintain Parallel Brace Sides

1. Slide the foot-bars together or apart as necessary to correspond with the narrowing or widening of thigh and calf cuffs. To avoid straining the knee or ankle joints, ensure that all side-bars remain parallel and vertical when viewed in the coronal plane, and are not being bent or twisted.
2. Tighten the foot-bar width adjustment screws to 5 Nm torque.

#### 4.4.5 Set the Device Knee Center Height

The device knee center height is adjusted by extending or retracting the lower extendable knee side-bars.

With the brace supported (e.g. on a table or on the floor) and Tectus in **Lock** mode, and the brace knee is fully extended:

1. Loosen, on the medial and lateral sides, the lower extendable knee side-bar's extension lock screws.

Note... Side-bar scales are read at the arrow position indicated where the row three dots on the fixed bar rail lines up with the scale on the extendable bar, see Fig. 15. The example shown is a medial lower side bar set at 3 cm above minimum extension. This equates to a device knee center height of:

42.5 cm (item B in Section 3.2, Fig. 3)  
 + 3 cm (side-bar reading)  
 = 45.5 cm above the foot sole , or  
 = 47.5 cm above the ground.

**It is critical that the lower extendable knee side-bars are set to the same height for the device to function correctly and to prevent damage to the knee and ankle joints.**


Fig. 15

2. Slide the medial and lateral extension bars at the same time out/in to increase/decrease respectively the brace knee height equally on the medial and lateral sides, until the knee and follower joint centers are at the same required vertical height (measurement C in Fig. 13 minus the calculated offset).
3. Check that the medial and lateral scale readings are equal (see Fig. 15); adjust as necessary, then tighten the extension lock screws to 5 Nm torque.

#### 4.4.6 Adjust Height of Upper Thigh Cuff and Strap (Optional Fit).

The upper thigh cuff and strap may be removed if not required, i.e. for shorter users with an inside leg measurement of less than ~62 cm, where using them would make the brace too high for their inside leg measurement at the medial cuff's lowest height setting.

Adjust the upper thigh cuff and strap heights with the device not worn, using the measured inside leg measurement and fine-adjust later when the device is worn.

1. If the upper cuff and thigh strap are required, attach as necessary (lateral upper cuff attaches via the lateral upper extendable side-bar, the medial side strap is free-floating).
2. Raise/lower the upper thigh cuff and strap by loosening the extension lock screws on the lateral upper extendable knee side-bar's, and slide the bar so that the upper thigh strap on the medial side is approximately 5 cm below the groin (equal to the user's inside leg measurement D in Fig. 13).
3. Tighten the upper side-bar extension lock screws to 5 Nm torque.

#### 4.4.7 Adjust Ankle Unit (Option)

1. If adjustment is available on the third part ankle unit, adjust the device ankle unit as required to set dorsiflexion/plantar flexion resistances (see third party ankle IFU).

#### 4.4.8 Adjust Brace Foot and Free Foot A/P Position and Foot Blocks

1. Loosen the heel stop sagittal position adjustment screw.
2. Move user's foot within brace foot anterior or posterior as necessary until the user's medial ankle center is as close as it can be to the device medial ankle center.
3. Slide the heel stop up to the user's shoe heel so they touch (or as close as possible), then tighten the heel stop screw to 5 Nm torque. Measure the new heel stop position relative to the fully extended position and adjust the heel stop on the free brace foot to match.
4. Moving the heel stop may affect the number of foot blocks required to support the user's foot. Slide in or out foot blocks from their open side to add or remove foot blocks to the slipper toe as necessary.
5. Replicate changes made to the foot block configuration on the free foot.

#### 4.4.9 Fit the Brace and Free Foot on the User and Check Knee Height

1. To aid fitting and adjustment, with Tectus knee module powered on and in default settings (see Section 4.4.1), select, as required, **Lock** mode or **Free** mode respectively to lock the brace knee or to allow the knee flex or extend easily.

## 2. Open all cuff and foot straps

Note...To open cuff straps, pull the cord on the magnetic CLICK fastener to lift and release the fastener, pull outwards on the BOA knobs to release the strap ratchet lock).

3. Loosen all the cuff A/P shift/rotation adjustment screws.
4. Loosen the heel stop adjustment screw and position the heel stop fully open.
5. Fit the brace around the user's affected leg, either while they are sitting with leg outstretched on the examination table, or seated on a chair.
6. Close and loosely fasten all straps.

Note: To tighten BOA™ straps, push in the BOA™ dial to engage the lock and rotate clockwise to tighten the strap.

7. Fit the device free foot to the shoe on the unaffected leg .
8. Set the device to **Lock** mode and ask the user to stand and hold on to support (e.g. parallel bars or examination table). The device will extend as they stand up but not flex in **Lock** mode.
9. Check the knee height is correct, adjust as necessary and retighten to 5 Nm torque.

### 4.4.10 A/P Align User's Knee Relative to Device Knee

A/P alignment of the user's knee position with respect to the device's knee A/P position is primarily achieved by adjusting the lower thigh cuff sagittal position and angle. However, shifting the thigh cuff position and moving the limb/knee accordingly will have a knock on effect on where the calf cuff sagittal position needs to be to allow correct limb positioning within the brace. Similarly, the heel stop position may need to shift.

1. With the device in **Lock** mode and the user either standing and supported, or sitting on the examination table, and with the affected leg in the brace fully extended, adjust (both sides equally) the lower thigh cuff sagittal A/P shift position and rotation (loosening straps as necessary). Aim to position the user's thigh such that the user's anatomical compromise knee center is at the correct A/P offset position relative to the device knee center (refer to Tectus Clinician's Manual 938465). Use laser to aid measurement/alignment.

Note at the same time as adjusting the lower thigh cuff A/P position, you may have to also adjust the A/P shift position and rotation of the calf cuff and further adjust the heel-stop position to allow correct overall positioning of the leg within the brace.

2. When the knee is correctly aligned, tighten the lower thigh cuff A/P position and rotation adjustment screws to 5 Nm torque.
3. Tighten the lower thigh cuffs straps by rotating the BOA® dials until secure, do not overtighten.
4. The upper side bars should remain near as possible vertical if standing, or if on parallel to the examination table surface.

### 4.4.11 Adjust Calf Cuff Sagittal Position and Angle of Rotation

1. Without moving the user's thigh further, adjust the calf cuff sagittal A/P position and angle of rotation until the cuff is touching and supporting the user's calf in the required position. If possible without shifting the knee A/P position alignment, set, as may be appropriate, up to a maximum of 10 degrees of anterior tilt of the knee relative to the ankle. Adjust the heel stop position if necessary at the same time. Replicate any heel stop repositioning made on the brace foot on the free foot heel stop.

2. When correctly aligned below the knee, tighten all cuff and heel stop adjustment screws to 5 Nm torque.
3. Fasten and tighten all the calf/shin cuff fasteners and BOA® straps and fasten and tighten all shoe straps to secure the brace in position. Do not overtighten.
4. If after tightening the straps the calf is not secure, loosen the straps and add packing pads as required between the user's clothed calf and the calf carbon cushion as necessary, then retighten the straps.

#### 4.4.12 Adjust Upper Thigh Cuffs Sagittal Position and Angle (Optional Fit)

1. If the upper thigh cuff is required and fitted to the brace, loosen the single upper thigh cuff adjustment screw to allow the thigh cuff and strap to slide horizontally and/or rotate as required in the sagittal plane.
2. Wrap the cuff strap around the users thigh, through the buckle and fold back to fasten with hook and loop fixings so that the strap and upper thigh cuff fit the user comfortably and the A/P position of the upper thigh adjustment screw is approximately vertically aligned with the lower thigh cuff adjustment screws A/P positions. Readjust the upper thigh cuff/strap height as necessary (see Section 4.4.6) and ensure there are no twists in the strap.
3. Tighten the adjustment screw to 5 Nm torque.

#### 4.4.13 Final Checks and Adjustments

1. Set the device to **Lock** mode If not already in **Lock** mode, and ask the user to stand if not already standing, and to hold on to parallel bars, examination table or chair for support. You must ensure they are given support as necessary until they have safe hold, are standing upright and the device is locked in full extension.
2. With user standing upright, check the knee A/P alignment with a vertical laser line, check the knee height alignment and ensure the device is orientated as required. Make any readjustments as necessary in the standing position.
3. After any readjustments, retighten all adjustment screws to 5 Nm torque, fasten and retighten all straps, ensuring the trouser leg or other clothing worn over the affected leg is folded neatly inside the brace to avoid discomfort.

**To prevent damage to clothing take care not to trap clothing between the moving parts of the device.**

The device is now fitted to the user and should be correctly aligned.

4. Set the device to **Sit** mode and ask the user to sit down on the chair carefully (holding onto chair or parallel bars for support). Note there may be too little or too much device resistance at this point as the Tectus has not yet been calibrated for the user, and may need manual adjustment via the app so the user can actually sit down.
5. Set the device to **Free** mode and gently flex and extend the knee and confirm that the device is not overly uncomfortable throughout the range of movement. If not, check the alignment again and readjust the device as necessary.
6. Set the device to **Lock** mode and ask the user to hold onto a chair or parallel bars for support and stand up. Help the user to stand as necessary.
7. The mechanical alignment is complete, now proceed to Tectus Knee Module Calibration and Configuration (see Section 5) to set up the Tectus knee module for the user.

---

## 5 Tectus Knee Module Calibration and Configuration

The Tectus Knee Module must be calibrated and configured specifically for each new user.

Please refer to the Tectus Clinician's Manual 938465 on how to calibrate and configure the Tectus Knee Module for a specific user.

**Advise the user to be aware at all times of the bulkiness of the brace and free foot so that they do not trip over their own feet or bump into nearby objects.**

Always check which knee extension assist spring is in the device before starting calibration, as the stiffness of that spring may not be appropriate for the new user's weight and activity level—swap out as necessary.

---

## 6 Device Operation

For device operation and all other Tectus-related information, including but not limited to device function, controls, operation, Tectus software application, battery charging and key fob programming information, please refer to the Tectus Clinician's Manual 938465.

For the required evaluation exercises, refer to the training notes.

On completion of the evaluation exercises, with the user supported, undo all straps and carefully remove the device. Reset the device using the app (See Section 7). If no other user is to be evaluated on the same day with this device, pack away the brace, brace foot and associated equipment in the Peli™ case ready for the next practitioner or user. If all evaluations are complete or the loan period is over, return the device in its Peli™ case it to Blatchford/your supplier.

---

## 7 Resetting the Device for the Next User

There is no specific reset feature on the app for the Clinician. Normally, after evaluating a user, it is intended that the device be returned to Blatchford/your supplier where it would be checked, reset and cleaned before being loaned out again.

However if you intend to evaluate more than one user with the device before returning the device to Blatchford, be aware that previous user's configuration data remains in the Tectus Knee Module until overwritten by resetting and reconfiguring the device using the app.

Therefore, after finishing evaluating a user, reset the following device parameters in the BASIC SETUP app screen as shown in Section 4.4.1 point 3, then power down the device so that it is ready for the next user.

Ensure that, for each new user, you follow the Device Set-up section of the Tectus Clinician's Manual 938465 and go through every calibration stage until fully calibrated for the new user so that all previous user data is overwritten.

If the default (as-supplied) knee extension assist spring (which is the stiffest spring) was swapped out for one of a different stiffness for the evaluation, always refit the default spring after the evaluation/before fitting the device to a new user.

Always clean/disinfect the brace and brace foot after use and before fitting to a new user, see Section 8 for details.

---

## 8 Maintenance

The device is not field-serviceable; it must be returned to Blatchford/your service provider either after each loan period or every 6 months, whichever period is shorter, or as otherwise defined in the loan agreement.

Maintenance must be carried out by competent personnel.

The device is on loan from Blatchford. If it develops any fault, becomes damaged or badly worn, stop using the device and contact Blatchford/your service provider for details.

It is recommended to clean the device and make the following checks before fitting to a new user:

- Check all construction screws are present and fastened tight.
- Check all adjustment screws are present.
- Check all extendable slide-bars are not damaged or distorted, and can slide up and down when extension lock screws are undone
- Ensure battery is fully charged
- Check Tectus Knee Module for
  - Signs of wear and tear
  - Hydraulic leakage
  - Security of the joint covers
  - No alarms present when switched on (such as Critical Service Indicator)

For any issues with the Tectus Knee Module, refer to the maintenance section of the Clinician's Manual 938465.

Note... The maintenance inspection form referenced in the Tectus Knee Module Clinician's Manual 938465 is not required for this device.

Any changes in performance of this device must be reported to Blatchford.

Changes in performance may include:

- Unexpected increase or decrease in knee stiffness
- Instability
- Reduced knee support (free movement)
- Any unrecognized beeps or warning lights
- Any unusual noise

### Cleaning and Disinfecting

The device must be disinfected using IPA (Isopropyl Alcohol) or Lysol after use/before fitting to a new user.

**Read and observe the safety warnings associated with the disinfectant used. Wear eye protection and protective gloves, and carry out cleaning in a well-ventilated area and/or use an air respirator.**

For general cleaning, use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers. Ensure water/liquids do not enter the charging socket. Thoroughly dry before use.

---

## 9 Limitations on Use

### **Intended Life**

A local risk assessment should be carried out based upon activity and usage.

### **Lifting Loads**

User weight and activity is governed by the stated limits.

Load carrying by the user is not allowed during evaluation.

### **Environment**

The device is intended for indoor use only. Avoid exposing the device to corrosive elements such as water, acids and other liquids except for cleaning and disinfecting materials (see Section 8). Also avoid abrasive contaminants such as soil and sand for example as these may promote premature wear.

Exclusively for use between  $-10^{\circ}\text{C}$  and  $50^{\circ}\text{C}$  ( $14^{\circ}\text{F}$  to  $122^{\circ}\text{F}$ ).


# 10 Assembly Instructions (Remove and Replace Cuffs)

1

T25

5 Nm


To remove the upper thigh cuff and strap, remove the single M5 x 16 mm Torx adjustment screw and remove the cuff and strap. Refit the screw in the cuff to prevent loss.

2

T25


5 Nm

Clamp Screw

To remove the lower thigh cuff, loosen the Tectus module clamp screw and rotate the Tectus module, then unscrew and remove four M5 x 16 mm Torx adjustment screws. Refit screws in the cuff to prevent loss.

3


Unplug battery power cable from the Tectus battery module, taking care not to damage pins. Remove the cuff including the battery module.

4

1.5

1.5 Nm


If replacement lower thigh cuff is not supplied with a Tectus battery module, remove battery module (remove three M2.5 x 20 mm screws) and fit to new cuff and tighten to 1.5 Nm torque.

5

T25

5 Nm

Remove calf cuff in same way as removing lower thigh cuff by unscrewing and removing four M5 x 16 mm adjustment screws.

6

Fit replacement cuffs in the reverse order of removal, noting that cuff widths may need to be adjusted to fit. Take care to fit screws through sliding mechanism(s) each side.

Lower thigh cuff only: After fitting the Tectus battery module to the cuff, reconnect the battery power cable to the battery module with cuff close to side-bars before fixing the cuff to the side bars, as it will be easier to access the connection port with the cuff free and help prevent potential damage to the connectors. See Clinician's Manual 938465 for correct connection orientation.

After connecting the cable, fix the calf cuff in position and secure hand tight with four M5 x 16 Torx adjustment screws.

# 11 Technical Data

## Device

Operating and Storage Temperature Range:	-10°C to 50°C (14°F to 122°F)
Charging Temperature Range:	0°C to 35°C (32°F to 95°F)
Device Weight (including Battery):	4.38 kg (9 lb 10.5 oz)
Activity Level:	(1), 2, 3, (4)
Maximum Combined User and Load Weight:	100 kg (220 lb)
Minimum User Weight:	50 kg (110 lb)
Brace Range of Motion at Knee:	110°
Cuff A/P Shift and Rotation Range:	
• Upper Thigh Cuff:	±25 mm, ±15°
• Lower Thigh Cuff:	±32.5 mm, ±7°
• Calf Cuff:	±32.5 mm, ±6°
Cuff Internal Width Ranges:	
• Standard Lower Thigh Cuff:	16.5–11.5 cm (6½–4½ in)
• Large Lower Thigh Cuff:	20.5–15.5 cm (8¼–6¼ in)
• Standard Calf Cuff:	15.6–10.6 cm (6¼–4¾ in)
• Large Calf Cuff:	18.6–14.6 cm (7¾–5¾ in)
Build Height:	See Section 3.2 <i>Height Adjustment Ranges</i>

## Device Battery

Battery Type:	Rechargeable Ni-MH 950 mAh 6V
Battery Life:	Variable, but max. 18 hrs from full charge. Charge daily.

## Battery Charger

Charging time to full charge:	2 hours
Input Voltage:	100V to 240V AC
Frequency:	50 Hz to 60 Hz
Max. Output Voltage:	21V DC
Max. Output Current:	0.35 A

## Key Fob Battery

Battery Type:	CR2032, Disposable Lithium 3V
---------------	----------------------------------

---

# 12 Ordering Information

## Order Example

<b>SU-TEB</b>	-	<b>L</b>
---------------	---	----------

Tectus  
Evaluation  
Brace

Left

<b>SU-TEB</b>	-	<b>R</b>
---------------	---	----------

Tectus  
Evaluation  
Brace

Right

## Transport & Storage

Ensure that the device is in a straight (fully extended) position and horizontal on a level surface when temporarily not in use, otherwise store it (flexed) in its Peli™ case.

If the device is to be stored or shipped, always store in its Peli™ case along with all tools and accessories as supplied by Blatchford.

If a unit is to be returned, contact Blatchford prior to sending if there is any sign of physical damage to the battery pack.

For transport purposes the unit must be switched off, and packed securely so that:

- a) It cannot be inadvertently switched on and become operative while in transit;
- b) The device is thoroughly protected from damage in transit.

# Manufacturer's Declarations and Legal Information

## Electromagnetic Environment

The device is designed for operation in the following electromagnetic environment:

- Operation in a professional healthcare facility (e.g. hospital) only.

## Electromagnetic Emissions

Interference Measurements	Compliance	Electromagnetic Environment Directive
HF emissions according to CISPR 11	Group 1/class B	The product uses HF energy exclusively for its internal functioning. Its HF emissions are therefore very low, and interference with neighboring electronic devices is unlikely.

## Electromagnetic Interference Immunity

Phenomenon	EMC Basic Standard or Test Procedure	Interference Immunity Test Level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
High-frequency electromagnetic fields	IEC 61000-4-3	10V/m 80 MHz to 6 GHz 80% AM at 400 Hz 100% AM at 1 kHz
	FCC CFR 47 part 15 (b) ANSI C63.4:2014	30 MHz to 12.5 GHz
	AIM Standard 7351731 (Rev 2.00) 02:2017 Radiated RF Immunity	433.92 MHz, 3 V/m 860-960 MHz, 54 V/m 2.45 GHz, 54 V/m
	AIM Standard 7351731 (Rev 2.00) 02:2017 Magnetic Field Immunity	134.2 kHz, 65 A/m 13.56 MHz, 12 A/m
Magnetic fields with rated power frequencies	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz

## Bluetooth® Module

IC Part No.	BGM220PC22HNA2
FCC ID	QOQ-GM220P
Frequency	2.4 GHz
Effective Radiated Power	8 dBm

---

## FCC Compliance Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Note... Contains pre-approved **Bluetooth**® module: FCC ID: QOO-GM220P IC: BGM220PC22HNA2.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation. This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party: Blatchford Inc.

1031 Byers Road, Miamisburg, OH 45342  
800-548-3534 (toll free) | 937-291-3636  
customerservice@blatchfordus.com

## RSS Compliance Statement

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website Blatchford Inc. 1031 Byers Road, Miamisburg, OH 45342.

800-548-3534 (toll free) | 937-291-3636 customerservice@blatchfordus.com

---

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## Warranty

This device is warranted for the duration of the loan period only.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to Blatchford.

## Environmental Aspects

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

## Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgments

Tectus and Blatchford are registered trademarks of Blatchford Products Limited.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG Inc., and any use of such marks by Blatchford is under License. Other trademarks and trade names are those of their respective owners.

The MRC Muscle Scale is used with the permission of the Medical Research Council, © Crown Copyright.

QuickFit, CLICK and BOA are registered trademarks of Click Medical 205 Hilltop Parkway, Warehouse 101, Steamboat Springs, CO 80487 USA.

Peli is a registered trademark of Peli Products S.L.U.

Velcro is a registered trademark of Velcro IP Holdings LLC

## Manufacturer's Registered Address

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Die MRC-Muskel-Skala wird mit Genehmigung des Medical Research Council verwendet, © Crown Copyright.
MRC Muscle-skalaen brukes med tillatelse fra Medical Research Council, © Crown Copyright.
**Blatchford Europe GmbH**