# Child's M30 4-Prong Laminate Adapter

## Instructions for Use

189626

EN	Instructions for Use	2
DA	Brugsanvisning	10
NO	Bruksanvisning	18
FI	Käyttöohjeet	26
SV	Bruksanvisning	34
RU	Инструкция протезиста	42
ZH	使用说明	50
AR	إرشادات الاستخدام	58

# Contents

EN

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	5
5 Maintenance .....	6
6 Limitations on Use .....	6
7 Fitting Advice.....	7
8 Socket Lamination Procedure.....	7
9 Technical Data .....	8
10 Ordering Information .....	8

# 1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout to refer to the Child's M30 4-Prong Laminate Adapter.

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively as part of a lower limb prostheses, intended for a single user.

This device is laminated into a prosthetic socket and connects via a Blatchford M30 Male Adapter, allowing rotational adjustment.

## Features

- Lightweight, strong, stainless steel construction

## Activity Level

This device is suitable for Activity Levels 1-4 (weight limits apply, see *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

Not recommended for short residua or where angular alignment changes might be required post-manufacture.

## Clinical Benefit

Provides an in-built connection point from the socket to other prosthetic componentry.

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted or excessive movement, or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, competitive running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling or running is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The user should be advised to contact their practitioner if their condition changes.


The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

## 3 Construction

### Principal Parts

- |                          |                 |
|--------------------------|-----------------|
| • Adapter                | Stainless steel |
| • Adapter Plug           | Nylon           |
| • M4x16 Socket Cap Screw | Plated steel    |
| • Lock Washer            | Plated steel    |

### Component Identification

## 4 Function

The rotation offered by the adapter allows for a customized limb alignment to suit the user. When the socket cap screw is tightened securely using Loctite 243 and the correct torque setting, this device offers a secure way to connect two parts of the limb.

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness, restricted/excessive movement or significant wear.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the socket cap screw for tightness, see *Construction* section. If loose, remove and clean, then apply Loctite and tighten to the correct torque setting
- Check for defects that could affect proper function

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear).

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

This device is waterproof to a maximum depth of 1 meter. Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water.

Exclusively for use between -15 °C and 50 °C  
(5 °F to 122 °F).

Suitable for submersion

## 7 Fitting Advice

*The instructions in this section are for practitioner use only.*

Symptom	Solution
A recurring noise occurs at the interface.	The Cap Screw must be tightened. Apply Loctite and tighten to the correct torque setting. Check for signs of lamination failure.
The adapter moves out of position.	Advise the user they must not use the device until adjusted, repaired or replaced. Inspect device. If components are worn or damaged, then remove and replace with new components. Reposition as necessary and retighten as per <i>Construction</i> section.

## 8 Socket Lamination Procedure

*The instructions in this section are for practitioner use only.*

Prior to lamination:

- Prior to shaping the prongs to suit the socket shape, turn back (unscrew) the 4-Prong Adapter a quarter of a turn to allow for adjustment of rotational alignment in the laminated socket
- Take care when shaping the prongs to suit the socket shape
- When shaping the prongs, ensure that the Blatchford M30 Male Adapter is fully engaged and that the Cap Screw is tightened


The lay up should be appropriate for the weight and activity of the user of the device.

Note... Remove the Adapter Plug after socket lamination.

## 9 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight:	100 g (3.5 oz)
Activity Level:	1–4
Maximum User Weight:	60 kg (132 lb)
Attachment Type:	M30 Child's Adapter
Range of Adjustment:	360° rotation
Build Height:	10 mm

### Build Height

### Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

## 10 Ordering Information

Device	Part Number
Child's M30 4-Prong Laminate Adapter	189626

## **Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## **CE Conformity**

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## **Compatibility**

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## **Warranty**

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

## **Reporting of Serious Incidents**

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## **Environmental Aspects**

This product is made from a recyclable material. Where possible, the components should be recycled in accordance with local waste handling regulations.

## **Retaining the Packaging Label**

The practitioner is advised to keep the packaging label as a record of the device supplied.

## **Trademark Acknowledgements**

Blatchford is a registered trademark of Blatchford Products Limited.

## **Manufacturer's Registered Address**

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**Endolite India Ltd.**
**Blatchford Inc.**
Blatchford Europe GmbH