

Epirus

Epirus incorporates an elastic Multiflex style ankle joint with an energy-storing-and-return prosthetic foot, which uses e-carbon foot springs to efficiently absorb energy during weight bearing and return it during off-loading, in order to aid propulsion. The elastic ankle joint provides a greater range-of-motion for improved adaptability. By providing plantarflexion motion at heel strike, cushioning/comfort improves compared to a fixed ankle ESR foot. The split-toe design permits further medial-lateral slope compliance.

Clinical Outcomes using e-carbon feet

Much research confirms the substantial equivalency of all energy-storing and return feet, including Blatchford e-carbon feet¹.

With respect to **SAFETY**

- High mean radius of curvature for Esprit-style e-carbon feet²: “The larger the radius of curvature, the more stable is the foot”

With respect to **MOBILITY**

- Allow variable running speeds³
- Increased self-selected walking speed⁴
- Elite-style e-carbon feet (L code VL5987) or VT units demonstrate the second highest mobility levels, behind only microprocessor feet⁵

With respect to **LOADING SYMMETRY**

- Users demonstrate confidence in prosthetic loading during high activity⁶
- Improved prosthetic push-off work compared to SACH feet⁷
- Increased prosthetic positive work done⁴

With respect to **USER SATISFACTION**

- High degree of user satisfaction, particularly with high activity users⁸

Clinical Outcomes using Multiflex-style ankles

Multiflex was the “habitual” foot for all or majority of participants in 13 different studies⁹⁻²¹

With respect to **SAFETY**

- Low stiffness at weight acceptance leads to early foot-flat and greater stability for lower mobility patients²²
- No loss of stability during standing with Multiflex than fixed ankle/foot²³
- Easier to walk on uneven ground with Multiflex than fixed ankle/foot^{23,24}
- Easier to walk up a slope with Multiflex than fixed ankle/foot²³

With respect to **MOBILITY**

- Little to no difference in gait mechanics compared to flexible, “energy storing” prosthetic feet²⁵
- Increased prosthetic ankle range-of-motion with Multiflex compared to fixed ankle/foot^{23,24,26-28}
- Increased prosthetic ankle power with Multiflex compared to fixed ankle/foot²⁴
- Prosthetic rollover shape closer to natural biomechanics than fixed ankle/foot²⁶
- Users can walk longer distances and report “smoother” gait with Multiflex compared to fixed ankle/foot²⁴
- Benefits in mobility for bilateral users^{23,24,26,27}
- Mixed objective results when user group was more active than is recommended for Multiflex^{29,30} so may benefit more from a similar but higher activity foot like Epirus.

With respect to **RESIDUAL LIMB HEALTH**

- Equivalent socket comfort to higher technology, energy-storing feet²⁹

With respect to **LOADING SYMMETRY**

- Improved stance phase timing symmetry with Multiflex compared to fixed ankle/foot²⁸
- Reduced sound limb loading with Multiflex compared to fixed ankle/foot²⁸

With respect to **USER SATISFACTION**

- Mixed subjective feedback around preferences when user group was more active than is recommended for Multiflex²⁵ so may benefit more from a similar but higher activity foot like Epirus.
- Majority of users rate Multiflex as either “good” or “acceptable”³¹ and prefer Multiflex to fixed ankle/foot²⁴

References

1. Crimin A, McGarry A, Harris EJ, et al. The effect that energy storage and return feet have on the propulsion of the body: A pilot study. *Proc Inst Mech Eng [H]* 2014; 228: 908–915.
2. Curtze C, Hof AL, van Keeken HG, et al. Comparative roll-over analysis of prosthetic feet. *J Biomech* 2009; 42: 1746–1753.
3. Strike SC, Arcone D, Orendurff M. Running at submaximal speeds, the role of the intact and prosthetic limbs for trans-tibial amputees. *Gait Posture* 2018; 62: 327–332.
4. Ray SF, Wurdeman SR, Takahashi KZ. Prosthetic energy return during walking increases after 3 weeks of adaptation to a new device. *J Neuroengineering Rehabil* 2018; 15: 6.
5. Wurdeman SR, Stevens PM, Campbell JH. Mobility analysis of Amputees (MAAT 5): Impact of five common prosthetic ankle-foot categories for individuals with diabetic/dysvascular amputation. *J Rehabil Assist Technol Eng* 2019; 6: 2055668318820784.
6. Haber CK, Ritchie LJ, Strike SC. Dynamic elastic response prostheses alter approach angles and ground reaction forces but not leg stiffness during a start-stop task. *Hum Mov Sci* 2018; 58: 337–346.
7. Rock CG, Wurdeman SR, Stergiou N, Takahashi KZ. Stride-to-stride fluctuations in transtibial amputees are not affected by changes in push-off mechanics from using different prostheses. *PLoS one*. 2018;13(10).

8. Highsmith MJ, Kahle JT, Miro RM, et al. Differences in Military Obstacle Course Performance Between Three Energy-Storing and Shock-Adapting Prosthetic Feet in High-Functioning Transtibial Amputees: A Double-Blind, Randomized Control Trial. *Mil Med* 2016; 181: 45–54.
9. Moore R. Patient Evaluation of a Novel Prosthetic Foot with Hydraulic Ankle Aimed at Persons with Amputation with Lower Activity Levels. *JPO J Prosthet Orthot* 2017; 29: 44–47.
10. Moore R. Effect on Stance Phase Timing Asymmetry in Individuals with Amputation Using Hydraulic Ankle Units. *JPO J Prosthet Orthot* 2016; 28: 44–48.
11. Buckley JG, De Asha AR, Johnson L, et al. Understanding adaptive gait in lower-limb amputees: insights from multivariate analyses. *J Neuroengineering Rehabil* 2013; 10: 98.
12. Sedki I, Moore R. Patient evaluation of the Echelon foot using the Seattle Prosthesis Evaluation Questionnaire. *Prosthet Orthot Int* 2013; 37: 250–254.
13. Rogers JP, Strike SC, Wallace ES. The effect of prosthetic torsional stiffness on the golf swing kinematics of a left and a right-sided trans-tibial amputee. *Prosthet Orthot Int* 2004; 28: 121–131.
14. Kobayashi T, Orendurff MS, Boone DA. Dynamic alignment of transtibial prostheses through visualization of socket reaction moments. *Prosthet Orthot Int* 2015; 39: 512–516.
15. Wright D, Marks L, Payne R. A comparative study of the physiological costs of walking in ten bilateral amputees. *Prosthet Orthot Int* 2008; 32: 57–67.
16. Vanicek N, Strike SC, Polman R. Kinematic differences exist between transtibial amputee fallers and non-fallers during downwards step transitioning. *Prosthet Orthot Int* 2015; 39: 322–332.
17. Barnett C, Vanicek N, Polman R, et al. Kinematic gait adaptations in unilateral transtibial amputees during rehabilitation. *Prosthet Orthot Int* 2009; 33: 135–147.
18. Emmelot C, Spauwen P, Hol W, et al. Case study: Trans-tibial amputation for reflex sympathetic dystrophy: Postoperative management. *Prosthet Orthot Int* 2000; 24: 79–82.
19. Boonstra A, Fidler V, Eisma W. Walking speed of normal subjects and amputees: aspects of validity of gait analysis. *Prosthet Orthot Int* 1993; 17: 78–82.
20. Datta D, Harris I, Heller B, et al. Gait, cost and time implications for changing from PTB to ICeX® sockets. *Prosthet Orthot Int* 2004; 28: 115–120.
21. de Castro MP, Soares D, Mendes E, et al. Center of pressure analysis during gait of elderly adult transfemoral amputees. *JPO J Prosthet Orthot* 2013; 25: 68–75.
22. Major MJ, Scham J, Orendurff M. The effects of common footwear on stance-phase mechanical properties of the prosthetic foot-shoe system. *Prosthet Orthot Int* 2018; 42: 198–207.
23. McNealy LL, A. Gard S. Effect of prosthetic ankle units on the gait of persons with bilateral trans-femoral amputations. *Prosthet Orthot Int* 2008; 32: 111–126.
24. Su P-F, Gard SA, Lipschutz RD, et al. Gait characteristics of persons with bilateral transtibial amputations. *J Rehabil Res Dev* 2007; 44: 491–502.
25. Boonstra A, Fidler V, Spits G, et al. Comparison of gait using a Multiflex foot versus a Quantum foot in knee disarticulation amputees. *Prosthet Orthot Int* 1993; 17: 90–94.
26. Gard SA, Su P-F, Lipschutz RD, et al. The Effect of Prosthetic Ankle Units on Roll-Over Shape Characteristics During Walking in Persons with Bilateral Transtibial Amputations. *J Rehabil Res Dev* 2011; 48: 1037.
27. Major MJ, Stine RL, Gard SA. The effects of walking speed and prosthetic ankle adapters on upper extremity dynamics and stability-related parameters in bilateral transtibial

amputee gait. *Gait Posture* 2013; 38: 858–863.

28. Van der Linden M, Solomonidis S, Spence W, et al. A methodology for studying the effects of various types of prosthetic feet on the biomechanics of trans-femoral amputee gait. *J Biomech* 1999; 32: 877–889.

29. Graham LE, Datta D, Heller B, et al. A comparative study of conventional and energy-storing prosthetic feet in high-functioning transfemoral amputees. *Arch Phys Med Rehabil* 2007; 88: 801–806.

30. Graham LE, Datta D, Heller B, et al. A comparative study of oxygen consumption for conventional and energy-storing prosthetic feet in transfemoral amputees. *Clin Rehabil* 2008; 22: 896–901.

31. Mizuno N, Aoyama T, Nakajima A, et al. Functional evaluation by gait analysis of various ankle-foot assemblies used by below-knee amputees. *Prosthet Orthot Int* 1992; 16: 174–182.