

Linx Amputee Profile Guide

A successful candidate for the Linx should:

- be able to fully load the prosthesis
- be able to walk without aid

Endolite PG# _____

Patient ID/Approval Number: _____

Age: _____ Sex: M F Side: L R

Anticipated Fitting Date: _____

Precise Body Weight: _____ lbs.
(without prosthesis, Max user weight 125kg (275 lb))

Measurements

Height: _____ ft. _____ in.

Residual Limb Length: _____ mm
Linx Build Height (w/o proximal)
470-565mm standard pylon (18 1/2" to 22 1/4")
up to 679mm extended pylon (26 3/4")

Present Prosthesis _____
(specify knee & foot)

Has true variable cadence: Y N

Able to negotiate down ramps without aid: Y N

Able to negotiate stairs, step over step: Y N


Insurance approval obtained: Y N

Bilateral: Y N

Occupation: _____

Activities: _____

Experience: _____



Foot Size: _____ cm Side: L R

Foot Shell: Light Dark

Facility Information

Prosthetist _____

Facility Name _____

Address _____

City _____

State/Zip _____

Phone _____

Fax _____

Email Address _____

Proximal Attachments

(Ordered separately, select choice)

Rotating Pyramid with Shift - 239017



Non-rotating Pyramid with Shift - 239089



Female Pyramid with Shift - 189128



M36 Threaded Adapter with Shift - 239092



Complete and return to determine if your patient is a successful candidate for the Linx.
fax: 800.929.3636 Attn: Customer Service or email: CustomerService@endolite.com

LINUXS042017-3FORM



endolite